



ASV Quick Guide



HAMILTON
MEDICAL
Intelligent Ventilation since 1983


 **Warning**

This Quick Guide is based on evaluations by clinicians within and outside Hamilton Medical and is intended to serve as an example. This Quick Guide does not replace either the official operator's manual of your ventilator or the clinical judgment of a physician. This Quick Guide should not - on its own - be used for clinical decision making.

Table of contents

1. Adaptive Support Ventilation (ASV) basics.....	4
2. Preparing for ventilation with ASV.....	6
3. Settings.....	8
4. ASV graph.....	12
5. Working principles of ASV.....	14
6. Monitoring ASV.....	16
7. Adjusting ASV.....	18
8. Weaning in ASV.....	20
Appendix.....	22
Lung-protective rules.....	22
Understanding the ASV safety box.....	24
Glossary.....	26

1. Adaptive Support Ventilation (ASV) basics



ASV focuses on simplification of mechanical ventilation. This means:

- ① Eliminating separate modes for passive and active patients
- ② Reducing controls relevant for CO₂ elimination to just %MinVol
- ③ Adding direct access to controls relevant for oxygenation (PEEP/CPAP and Oxygen)

ASV maintains an operator-preset minimum minute volume independent of the patient's activity. The breathing pattern (tidal volume, rate, and inspiratory time) is calculated, based on the assumption that the optimal breathing pattern results in:

- a The least work of breathing
- b The least amount of ventilator-applied inspiratory pressure

A **lung-protection strategy** ensures ASV's safety. ASV attempts to guide the patient using a favorable breathing pattern and avoids potentially detrimental patterns such as rapid shallow breathing, excessive dead space ventilation, intrinsic PEEP, barotrauma, and volutrauma.

2. Preparing for ventilation with ASV

The image shows a ventilator control panel interface. At the top right, the mode is set to **ASV** (Adult/Ped.). The main display area includes a **Standby** button (highlighted in yellow), a patient selection dropdown set to **Adult/Ped.**, and a **Last patient** button. Below these are buttons for **ASV**, **NIV**, and **CPR**. Further down are **Male** and **Female** gender selection buttons, a patient height of **174 cm**, and a weight of **70^{IBW} kg**. At the bottom of the main area are **Preop check** and **Start ventilation** buttons. A note at the bottom of this section reads: "To power off, press power/standby key > 3 s". On the right side, there are three circular gauges: **100%** for **%MinVol**, **5 cmH2O** for **PEEP/CPAP**, and **50%** for **Oxygen**. Below these gauges are **Controls** and **Alarms** buttons. At the bottom of the screen are **Monitoring**, **Utilities**, **Events**, and **System** buttons, along with a power status indicator showing **INT** and **AC**.

Preparation

- Perform pre-operational check of the ventilator.
- Select accurate patient height and gender for the calculation of the ideal body weight (IBW).

3. Settings

ASV initial setup

ASV Adult/Ped.

Modes: Standby (highlighted), Adult/Ped., Last patient

1 100 % \dot{V}_{MinVol}

2 5 cmH₂O PEEP/CPAP

3 50 % Oxygen

Male Female 174 cm Pat. height 70 kg IBW

5 Start ventilation 4 Controls

To power off, press power/standby key > 3 s

Alarms

Monitoring Utilities Events System

Basic More TRC

30 cmH₂O Pasp/limit

174 cm Pat. height

IBW: 70 kg MinVol: 7.0 l/min

5.0 l/min Flowtrigger

1 100 % %MinVol

2 5 cmH₂O PEEP/CPAP

3 50 % Oxygen

4 Controls

Alarms

Basic More TRC

50 s P-ramp

25 % ETS

Sigh

4 Controls

Alarms

5 Start ventilation

- 1 **%MinVol:** Suggested initial setting for a normal patient: 100% (ARDS: 120%)
For adults, minute volume is calculated at 0.1 l per kg of IBW. For a patient with IBW = 70 kg, 100% MinVol results in 7 l/min, 50% MinVol is 3.5 l/min, 200% MinVol is 14 l/min.
For pediatric patients, minute volume is calculated in a range from 0.3 l per kg for IBW = 3 kg to 0.1 l per kg for IBW = 30 kg.
- 2 **PEEP/CPAP:** Suggested initial setting: 5 cmH₂O (or according to your ICU standard)
- 3 **Oxygen:** Suggested initial setting: 50% (or according to your ICU standard)
- 4 **Controls:** In the Controls window, check the default settings. If required, adjust the following settings according to the patient's condition:
 - Maximum pressure set by ASV (**Pasvlimit**). Default for normal patients: 30 cmH₂O
 - Flow or pressure trigger
 - Pressure ramp (**P-ramp**)
 - Expiratory trigger sensitivity (**ETS**)
- 5 Connect the patient to the ventilator and touch **Start ventilation** to start

3. Settings


Alarms

Limits 1

Parameter	Value
Pressure	40
ExpMinVol	10
fTotal	40
Vt	850
Apnea time	20

PEEP/CPAP: 5 cmH₂O
Oxygen: 50%

Buttons: Auto, Controls, Alarms



Modes ASV Adult/Ped.

Parameter	Value
Basic	30 cmH ₂ O
Pasvlimit	174 cm
Pat. height	5.0 l/min
%MinVol	100%
PEEP/CPAP	5 cmH ₂ O
Oxygen	50%

Flowtrigger: 5.0 l/min
IBW: 70 kg
MinVol: 7.0 l/min


Buttons: More, TRC, Controls, Alarms

Check that the high Pressure alarm limit is set to an appropriate value.

Suggestions: Normal patient: 40 cmH₂O for 100% MinVol


- i** The maximum inspiratory pressure delivered in ASV (**Pasvlimit**) will be **10 cmH₂O below the preset high pressure limit**, indicated by a blue band on the pressure curve graph. The maximum inspiratory pressure for ASV can be also set using the **Pasvlimit** control in the Controls window. Changing the **Pasvlimit** value also changes the high Pressure limit.
- !** To avoid lung over distension, **check the Vt high alarm limit**, and make sure the target minute ventilation can still be reached in **passive** patients. Inspiration is aborted in mechanical breaths as soon as the volume exceeds 1.5 x Vt high alarm limit.

4. ASV graph




- ① Horizontal axis for respiratory rate (f)
 - ② Vertical axis for tidal volume (V_t)
 - ③ Minute volume curve – see next page
 - ④ Safety frame in which target point may move
 - ⑤ Target point, formed by intersection of target tidal volume and target rate
 - ⑥ Current patient values, formed by intersection of measured tidal volume and current respiratory rate
-
- ① ASV adjusts the settings to guide the patient's current values to the target point. This may be achieved for spontaneous breathing or passive patients. When the patient condition meets the target, the patient is considered optimally ventilated according to ASV. **It is not an indication of the patient's clinical condition.**

5. Working principles of ASV



Optimal combination of tidal volume / respiratory rate
(in this example 15 x 400 ml for a minute volume of 6l)

Taking all possible combinations of respiratory rate and tidal volume into account, ASV calculates the optimal breathing pattern based on operator entries of %MinVol and the IBW, as well as on the measurement of RC_{exp} . The device works on the assumption that the optimal breath pattern is identical to the one a totally unsupported patient will choose naturally (least work of breathing).



Lung-protective rules

ASV applies a lung-protective rules strategy to avoid

- Ⓐ High tidal volumes and pressures
- Ⓑ Low alveolar ventilation
- Ⓒ Dynamic hyperinflation or breath stacking
- Ⓓ Apnea

This lung-protection strategy ensures ASV's safety while it maintains an operator-preset, minimum minute ventilation independent of the patient's activity.

See Appendix I for detailed rules.

6. Monitoring ASV

Expiratory time constant (RC_{exp})

The expiratory time constant (RC_{exp}) is a measure of how fast or slow the lung fills and empties. It is the product of compliance and resistance. Thus, this simple measurement assesses the two main characteristics of respiratory mechanics.

Why monitor RC_{exp} ?


- RC_{exp} is used as an input in ASV: if ASV selects an unexpected tidal volume-respiratory rate combination, checking RC_{exp} helps to understand why.
- To understand the respiratory mechanics of the patient.
- To set the breath cycle: To have a complete expiration, expiratory time should be at least equal to $2 \times RC_{exp}$

How is RC_{exp} measured?


RC_{exp} is measured breath-by-breath as the ratio between volume and flow during expiration. It is displayed in the Monitoring window, and is accurate in all breaths with passive exhalation.

Normal values in intubated adult patients


- 1 Short: < 0.6 s: restrictive disease: ARDS, atelectasis, chest wall stiffness
- 2 Normal: $0.6 - 0.9$ s: normal compliance and resistance or combination of decreased compliance and increased resistance
- 3 Long: > 0.9 s: obstructive disease (COPD, asthma), bronchospasm, endotracheal tube obstruction or incorrect positioning



1



2



3

7. Adjusting ASV

Management of %MinVol

Passive patient

- If PaCO₂ is too high, increase %MinVol
- If PaCO₂ is too low, reduce %MinVol




Active patient

- If the patient is tachypneic and/or has a high respiratory effort, increase %MinVol
- If the patient's respiratory rate is lower or pressure support (Pinsp) is higher than desired, decrease %MinVol

- i** Check blood gas analysis after 30 min, and adjust as needed
- i** The required %MinVol may be as high as 200%, but should not exceed 300%. If a very high %MinVol is required to meet the ventilation demand, consider using or increasing sedation, and control other factors such as high fever, sepsis, or metabolic acidosis. Except for some special cases such as hypothermia or chronic hypercapnea, a %MinVol lower than 100% is appropriate only for promoting spontaneous breathing activity.

Patient condition	%MinVol Adjustment
Patient is passive and has a high PaCO₂ and very low arterial pH .	Increase %MinVol in steps of 10% – 20% to restore normal PaCO ₂ and arterial pH.
Patient is active and continuously shows signs of respiratory distress .	Increase %MinVol in steps of 20%. After every step, wait for 2 to 5 minutes to observe the patient's response.
Patient is active , but spontaneous breathing disappears or becomes irregular .	Reduce %MinVol in steps of 10%. Observe patient's response after every adjustment.

8. Weaning in ASV



Oxygenation		CO2 elimination		Spont/Activity	
40	8	10.5	10	100	75
21	0	3.5	0	10	100
00:01	00:01	00:00			
Oxygen	PEEP	MinVol	Pinsp	RSB	%Spont
30	5	7.6	26	---	0
%	cmH2O	l/min	cmH2O	1/(l*min)	%

Oxygenation		CO2 elimination		Spont/Activity	
40	8	10.5	10	100	75
21	0	3.5	0	10	100
00:08	00:08	00:06	00:00	00:01	00:02
Oxygen	PEEP	MinVol	Pinsp	RSB	%Spont
21	5	8.7	10	26	100
%	cmH2O	l/min	cmH2O	1/(l*min)	%

Phase 1 - Screening

If deep sedation is stopped and the patient is active, gradually reduce %MinVol (at most to 70% MinVol), PEEP, and Oxygen every hour.

Phase 2 - Observation

If the patient's respiratory rate is < 30 breaths/min, $\text{Pinsp} < 15 \text{ cmH}_2\text{O}$, $\text{PEEP} \leq 8 \text{ cmH}_2\text{O}$, $\text{Oxygen} \leq 40\%$ or according to your ICU standard for 30 min to 2 h, consider an SBT.

Phase 3 - Spontaneous breathing trials (SBT)

Perform SBT. Suggested SBT settings:

PEEP = $5 \text{ cmH}_2\text{O}$, Oxygen = 30%

%MinVol = 25% for 30 minutes

If SBT is successful after 30 minutes and extubation criteria are fulfilled, consider extubation.

- i** This protocol is our suggestion. You may have your own specific weaning criteria that patients have to fulfill before starting SBTs.

Appendix I

Lung-protective rules

A) High tidal volume limit

The tidal volume applied by ASV is limited by three operator settings: P_{asvlimit} , V_t high alarm limit, and IBW. The maximum pressure to be applied in the ASV mode is 10 cmH₂O below the high Pressure limit. Maximum tidal volume is limited by $(P_{\text{asvlimit}} - PEEP) \times \text{compliance}$. In addition target volume is limited to $1.5 \times V_t$ high limit, and pressure support is limited in such a way that the inspired volume does not exceed the V_t high limit.

B) Low tidal volume

It is widely accepted that a first approximation of dead space can be obtained by the following simple equation (Radford 1954): $V_{Daw} = 2.2 \times \text{IBW}$. The lower limit for tidal volume is based on this equation and calculated to be at least twice the dead space. The minimum V_t is $4.4 \times \text{IBW}$.

C) High rate limit

The equation used to calculate the maximum rate is: $f_{\text{max}} = \text{target MinVol} / \text{minimum } V_t$.

To achieve a nearly complete exhalation to the equilibrium point of the respiratory system (90% of the maximum potential volume change), an expiratory time of at least $2 \times RC_{\text{exp}}$ is theoretically

required. For this reason, ASV calculates the maximum rate based on the principle of giving a minimum inspiratory time equal to $1 \times RC_{exp}$ and a minimum expiratory time equal to $2 \times RC_{exp}$, which results in these equations:

$$f_{max} = 60 / (3 \times RC_{exp}) = 20 / RC_{exp}$$

$$f_{max} \leq 60 \text{ b/min}$$


ASV always uses the lower of the two values. This limit applies to the respiratory rate of the ventilator only, not to the respiratory rate of the patient.

D) Low rate limit


The lowest target rate for adult patients is fixed at 5 b/min. For pediatric patients, the lowest target rate is in a range from 7.5 b/min for IBW = 30 kg to 15 b/min for IBW = 3 kg.

Appendix II


Understanding the ASV safety box




Wide square-shaped box: normal lung mechanics ($RC_{exp} = 0.6$ s).



Low and wide box: low compliance or a «stiff» lung ($RC_{exp} = 0.3$ s)



Narrow and high box: obstructions with long time constants and high resistance ($RC_{exp} = 1.2$ s)



Low and narrow box: high resistance and low compliance ($RC_{exp} = 0.8$ s)

Appendix III

Glossary

ASV	Adaptive Support Ventilation
ETS	Expiratory trigger sensitivity
FiO ₂	Fraction of inspired oxygen
%fSpont	Spontaneous breath percentage
PaCO ₂	Arterial partial pressure of CO ₂
Pasvlimit	Maximum pressure set by ASV, always equal to pressure limit -10.
PEEP	Positive end expiratory pressure

P _{insp}	Inspiratory pressure
%MinVol	Target minute volume on IBW ratio
RC _{exp}	Expiratory time constant
RC _{insp}	Inspiratory time constant
RR	Number of breaths per minute
SBT	Spontaneous breathing trial

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland

 +41 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

P/N 610684/03

Specifications are subject to change without notice.
© 2014 Hamilton Medical AG. All rights reserved.



HAMILTON
MEDICAL
Intelligent Ventilation since 1983