Instructions for use

PN 612375/04 2021-03-25

Nasal cannula for high flow oxygen therapy

REF	282495 Nasal cannula, size S, Adult/pediatric 282496 Nasal cannula, size M, Adult/pediatric 282497 Nasal cannula, size L, Adult
QTY	Quantity
	Follow the Instructions for use
	Do not reuse
	Applicable to adult patient group
	Applicable to pediatric/adult patient groups
	Do not use if the packaging is damaged
	Do not use any blades, knives, or cutters to open; they can damage the product
MD	Medical device

WARNING

- Reprocessing single-use products can affect physical and chemical material properties and may endanger patients. It can lead to increased particle release, an increased tendency to crack, and have adverse effects on microbiological resistance. Furthermore, an altered surface structure may result in a microbial aggregation of spores, allergens, and pyrogens, for example, or cause an increase in the number of particles released as a result of chemical changes in the material properties.
- Visually inspect the parts for damage. Discard if there is any sign of damage or any changes are identified.
- Used nasal cannulas may be contaminated. Observe the infection control procedures at your clinic. Follow applicable local laws and regulations regarding waste management.
- The manufacturer or distributor does not warrant the proper functioning of nasal cannula if they are cleaned or reused by the user.

CAUTION

- (USA only) Federal law restricts this device to sale by or on the order of a physician.
- This product is for use only by trained medical personnel or under the direction of qualified medical personnel.
- DO NOT use the nasal cannulas without humidification as it may lead to patient airway damage.
- DO NOT fully obstruct nostrils as it may lead to high pressure and patient injury.
- Keep the device away from flames and sparks.

NOTE

- Replace the nasal cannula in accordance with the applicable infection control policies and procedures, or if deemed necessary due to secretions, deposits, or nebulized drug therapy.
- Read the operating instructions for the appropriate high flow oxygen therapy device.
- We recommend using a water trap when using a collapsible tube, as condensation can occur when a heated wire is not used.
- The maximum delivered flow depends on the specifications and performance of the oxygen therapy device. In some markets, the maximum possible flow setting may be limited.

Intended use

The nasal cannula is a two-pronged device intended for direct flow administration of oxygen via the nasopharyngeal route.

This product is intended for single use only, for a maximum of 7 days.


Setting up the breathing circuit

To connect to the breathing circuit

- ▶ Connect the nasal cannula to the breathing circuit.

To fasten the nasal cannula


- 1 Put on a neck lanyard and nasal cannula.

- 2 Adjust the tightness of the nasal cannula.

- 3 Connect the nasal tube and breathing circuit.

- 4 Connect the tubing to the humidification chamber of the heated humidifier. Follow the instructions provided with the humidifier.

To secure the tubes

- ▶ Make sure the flow sensor measurement tubes are securely connected to the breathing circuit.

Technical specifications

Oxygen concentration	21% to 100%
Weight	size S, M, L: 40 g
Material	Silicone, PP, EVA, nylon All materials meet ISO 10993-5, <i>Biological evaluation of medical devices</i> .
Connection	OD22/ID15
Tube length	25 cm (9.8 in)
Operating time	Maximum 7 days
Regulatory	The manufacturer has established and applies a quality management system and fulfills ISO 13485. The nasal cannulas have been designed to comply with Council Directive 93/42/EEC.

Operation and storage conditions

Operating	Temperature:	-10°C to 50°C (14°F to 122°F)
Storage	Temperature:	-10°C to 50°C (14°F to 122°F)
	Humidity:	30% to 95%

Technical Data

PN	Size	Prong diameter (mm)	Inter-prong distance (mm)	Prong length (mm)
282495	S	4.0	11.8	13
282496	M	5.0	12.4	11
282497	L	6.0	11.1	13

PN	Size	Flow (l/min) at resistance of 2 mbar	Max. flow (l/min) at a back pressure of 50 mbar
282495	S	11.16	60
282496	M	16.19	80
282497	L	18.68	80

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